



**Argonne**  
NATIONAL  
LABORATORY

*... for a brighter future*



U.S. Department  
of Energy

UChicago ►  
Argonne<sub>LLC</sub>



**Office of  
Science**  
U.S. DEPARTMENT OF ENERGY


A U.S. Department of Energy laboratory  
managed by UChicago Argonne, LLC

# **MCP and Photocathode Testing and Systems Integration At the Advanced Photon Source**

*B. Adams, K. Attenkofer, M. Chollet, Z. Insepov, J. McPhate,  
O. Siegmund, D. Walters, M. Wetstein, Z. Yusof  
for the LAPPD Collaboration*




# Goals of the APS Test Stand






# Goals of the APS Test Stand





# Goals of the APS Test Stand






# The Current Setup

- Vacuum chamber operating at  $10^{-7}$  torr level
- MCP/photocathode assembly mounted on optical cage system.
- Cage system attached to side-mounted flange with SMA and HV feedthroughs.
- Operation with or without photocathode (CsI on diamond)
- Ti:Sapphire laser (50 fs, 800 nm), frequency-tripled to 266 nm
- Voltage on photocathode: 0 - 4.0kV
- Voltage on MCP from 1.5-2.0kV
- Timing measurements using 8-GHz and 16 GHz scopes


## Vacuum Cross





# The Current Setup



## Cage system with MCP









# The Current Setup







## Results: Photocathode Measurement









# Results: MCP Measurement





- Characterization of commercial Photonis MCP (The Chevron Model 3025 ).
- Amplification measured as integrated charge (# electrons) collected on a single stripline.
- Expected (total) amplification at 2kV:  
 $\sim 1 \times 10^7$




Example avg. scope signals at different MCP voltages.





## Results: MCP Measurement






# Results: MCP Measurement







## Results: MCP Measurement



- After characterizing the Photonis MCP, we coat the plates with 10 nm  $\text{Al}_2\text{O}_3$ .
- The “after-ALD” measurements have been taken without scrubbing.
- These measurements are ongoing.



# Coming Soon: Tests of Functionalized Borosilicate Samples





# Near Future Plans:

## The 'B' Configuration


- More compact arrangement of MCP's directly against a single flange. Minimal or no cabling. Simple or no photocathode. Simple stripline structure.
- Used for a precise and direct comparison of single or double channel plates, with all other variables held as constant as possible.
- Designed for simplicity, vacuum compatibility, interchangeability.
- Can be built while measurements are still taken on the current setup.
- Optical setup built onto modular, portable breadboards, and designed to handle a wide range of light sources.
- Can also be used with a well defined commercial MCP for photocathode characterization.
- Can be docked with a larger vacuum transfer system.







# Near Future Plans:

## The 'B' Configuration




(Can also be attached to a single flange with glass window to form a compact, two-flange "MCP")






# MCP Holder




Courtesy of Jason McPhate






# 6" Con-flat flange: (Side Feedthroughs)

Bottom




- 4 High Voltage Feedthroughs
- 1 Feedthrough for pumping down

Top




Courtesy of Dean Walters





# 6" Con-flat flange:

## 3 Basic Types

Flange with striplines for precision timing measurements



Flange with a single plate for amplification measurements



Flange with phosphor screen for testing image uniformity of single plates



# Summary

- We have successfully assembled the right resources, man-power, expertise, and experience necessary to meet our testing goals.
- We are presently following 2 parallel tracks:

## Current Setup

- Finishing up characterization of commercial MCP, before and after ALD.
- Long term gain study of MCP after ALD coating.
- Proof-of-principles test of MCP made using borosilicate glass with ALD coating.

## Future Setup

- Finishing up design phase.
- Plan to be ready for comprehensive testing of ALD-based channel plates within the next two months.
- Working on developing vacuum transfer capabilities for tests of photo-cathode samples.