



... for a brighter future

Fast X-Ray Detectors for the APS and Elsewhere

Bernhard W. Adams¹

¹Advanced Photon Source
Argonne National Laboratory

2010-06-11



UChicago ►
Argonne_{LLC}



A U.S. Department of Energy Laboratory
managed by UChicago Argonne, LLC

Contents

X-Ray Detector Needs

What is already out there

Synchrotron

Medical

Application examples

MCP or columnar structure as photocathode

A positively gated photocatode for nuclear resonance

Contents

X-Ray Detector Needs

What is already out there

Synchrotron

Medical

Application examples

MCP or columnar structure as photocathode

A positively gated photocathode for nuclear resonance

X-Ray Detector Needs

- ▶ Synchrotron: 2-d spatial and/or time resolution
- ▶ XFEL: 2-d spatial resolution, tons of x-rays
- ▶ medical: imaging, PET
- ▶ homeland security: large-area imaging
- ▶ some special applications

Contents

X-Ray Detector Needs

What is already out there

Synchrotron

Medical

Application examples

MCP or columnar structure as photocathode

A positively gated photocathode for nuclear resonance

What is Already Out There

- ▶ CCDs: direct-absorption and scintillator; can be gated, in principle, to μs or so, but not usually done. Spatial resolution down to a few microns, without gating down to ms resolution, \$\$ / cm^2
- ▶ Pixel-array detectors, can be gated to 100 ns, spatial resolution ca. 100 microns, without gating down to μs resolution, depending on ROI, \$\$ / cm^2
- ▶ wire chamber (gas avalanche amplification): location/time-tagged events μs resolution, area is cheap
- ▶ Our detector produces a **stream of location-time-tagged events, not just events in a gate interval**. Like wire chamber, but faster/higher rate and better x-ray stopping power, area is cheap

Contents

X-Ray Detector Needs

What is already out there

Synchrotron

Medical

Application examples

MCP or columnar structure as photocathode

A positively gated photocathode for nuclear resonance

Synchrotron

- ▶ Synchrotron: 2-d spatial and/or time resolution
- ▶ wide range of event rates, possibly spatially highly inhomogeneous photon load
- ▶ use scintillator with visible-light photocathode: slow response or low efficiency
- ▶ planar, direct-conversion photocathode: fast, low efficiency
- ▶ NB: at synchrotron, 100 ns time resolution is as good as 100 ps (bunch duration), faster than 100 ps very valuable

Contents

X-Ray Detector Needs

What is already out there

Synchrotron

Medical

Application examples

MCP or columnar structure as photocathode

A positively gated photocathode for nuclear resonance

Medical

- ▶ High detection efficiency at higher photon energies: 30 .. 150 keV, or even 511 keV for PET
- ▶ medical x-ray needs “no” time resolution, absorption-edge angiography could profit from ca. 1 ms, PET needs 100 ps or better
- ▶ large area

Contents

X-Ray Detector Needs

What is already out there

Synchrotron

Medical


Application examples

MCP or columnar structure as photocathode

A positively gated photocatode for nuclear resonance


Application Example: Protein Crystallography

white-beam Laue diffraction: one crystal, many λ



rotate crystal
to obtain series
of diffraction maps
intensities in diffraction s
FT into structure

LAPPD instead of CCD:
no readout overhead
rapid rot., continuous readout




Example: Dynamic Light Scattering of X-Rays


time evolution of spatial correlations

measure correlation of speckle patterns at 2 times


homodyne speckle



heterodyne speckle



challenge: wide dyn. range



M. Sutton et al., Opt. Expr. **11**, 2268-2277 (2003)


Example: XFEL: Tons of Photons in a Pulse

Problem at XFEL: Detector overload by huge pulse

solution: spread flux in time and space

in time: use efficient, slow scintillator

in space: grazing incidence



Contents

X-Ray Detector Needs

What is already out there

Synchrotron

Medical

Application examples

MCP or columnar structure as photocathode

A positively gated photocatode for nuclear resonance



MCP or columnar structure as photocathode

Leverage MCP/ALD technology developed in LAPPD

electron escape depth ca. 100 nm



x-rays mostly pass through CsI

10 keV \rightarrow ca. 1% abs.



use grazing incidence

\rightarrow MCP or nanocolumns
with ALD'd photocathode



Contents

X-Ray Detector Needs

What is already out there

Synchrotron

Medical

Application examples



MCP or columnar structure as photocathode

A positively gated photocatode for nuclear resonance

Nuclear Resonance, Time Domain


forward scattering from foils of increasing thickness

v. Bürck et al., PRB **46**, 6207-6211 (1992)



How this is measured on a synchrotron

Using synchrotron radiation: pulsed, broadband source




Gated MCP

new idea:


use gating pulse to actively clear out MCP
for 100 ps during SR pulse

MCP is ready 100 ps after SR pulse,
no “afterglow”



100-ps gated MCP

Katayama et al.,
Rev. Sci. Instrum. **62**, 124 (1990)



10^{-10} contrast possible?